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Ophthalmologiste Praticien Hospitalier
CHU Saint-Étienne – France

Mots clés

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


Perte entre 75 et 100%

Bibliographie

1. Parasympathetic Innervation of the Meibomian Glands in Rats - Mark S. LeDoux et al. - *Investigative Ophthalmology & Visual Science*, October 2001, Vol. 42, No. 11
2. Characterization of the innervation of the meibomian glands in humans, rats and mice - Bründl, M. et al. *Annals of Anatomy* (2021), Vol. 233.
3. Neurotransmitter Influence on Human Meibomian Gland Epithelial Cells - Wendy R. Kam and David A. Sullivan - *Investigative Ophthalmology & Visual Science*, November 2011, Vol. 52, No. 12
4. The Dopaminergic Neuronal System Regulates the Inflammatory Status of Mouse Lacrimal Glands in Dry Eye Disease - Ji, Yong Woo et al. *Investigative Ophthalmology & Visual Science* (2021), Vol. 62.
5. TFOS DEWS II Tear Film Report - Willcox Mark et al. - *The Ocular Surface*, 2017, p 441 to 515
6. The neurobiology of the meibomian glands - Cox SM, Nichols JJ - *Ocular Surface*, July 2014
7. Multicenter Study of Intense Pulsed Light Therapy for Patients with Refractory Meibomian Gland Dysfunction - Reiko Arita, et al. - *Cornea* Volume 37, Number 12, December 2018
8. Rosacea: Molecular Mechanisms and Management of a Chronic Cutaneous Inflammatory Condition - Yu Ri Woo, et al. *International Journal of Molecular Sciences*, September 2016
9. Rosacea: Epidemiology, pathogenesis, and treatment - Barbara M. Rainer et al. - *DERMATO-ENDOCRINOLOGY* 2018, VOL. 9, NO. 1, e1361574 (10 pages)
10. Treatment of ocular rosacea - Edward Wladis et al. - *Survey of Ophthalmology* (2018), Vol.63.
11. Improved telangiectasia and reduced recurrence rate of rosacea after treatment with 540 nm-wavelength intense pulsed light: A prospective randomized controlled trial with a 2-year follow-up - Luo, Y. et al. - *Experimental and Therapeutic Medicine* (2020), Vol. 19.
12. Therapeutic Effect of Intense Pulsed Light on Ocular Demodocosis - Zhang, X., et al. - *Current Eye Research* 2019, Vol. 3.
13. Intense Pulsed Light Therapy for Patients with Meibomian Gland Dysfunction and Ocular Demodex Infestation - Cheng et al. - *Current Medical Sciences* (2019), Vol.39.
14. Long-term effects of intense pulsed light treatment on the ocular surface in patients with rosacea-associated meibomian gland dysfunction - Seo Kyoung Yul et al. - *Contact Lens and Anterior Eye* (2018), Vol. 41.
15. TFOS DEWS II Tear Film Report - Willcox Mark et al. - *The Ocular Surface* (2017), Vol.15.
16. Intense Pulsed Light for the Treatment of Dry Eye Owing to Meibomian Gland Dysfunction - Vigo, L. et al. - *Journal of Visualized Experiment* (2019), N°146.
17. Meibum Expressibility Improvement as a Therapeutic Target of Intense Pulsed Light Treatment in Meibomian Gland Dysfunction and Its Association with Tear Inflammatory Cytokines - Choi, M. et al. - *Scientific Reports* (2019), Vol.9.
18. TFOS DEWS II pain and sensation report - Belmonte Carlos, et al. - *The Ocular Surface* (2017), Vol.15.
19. Analysis of Cytokine Levels in Tears and Clinical Correlations After Intense Pulsed Light Treating Meibomian Gland Dysfunction - LIU, R et al. - *American Journal of Ophthalmology* (2017).
20. Effect of inflammation on lacrimal gland function - Driss Zoukhri - *Experimental Eye Research*, May 2006; 82(5): 885-898
21. Aqueous deficiency is a contributor to evaporation-related dry eye disease - Charles W. McMonnies - *Eye and Vision* (2020) 7:6.

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